



✦ Extracting Forceps

✦ Root Elevators



01
SECTION



We will Provide You Extracting Forceps English & American Pattern in these Handles



FEATURE PRODUCTS

- A01: Manufactured from AISI 240 German stainless steel.
- A02: Highly polished finish for aesthetic and cosmetic resistance.
- A03: Smooth handle knurling for better grip and control.
- A04: Cross serrated beaks for strong and positive tooth grip.
- A05: Teflon washer in the joints for years of trouble free operation.
- A06: Fully autoclavavle/reusable.

Our delivery ready in forging stock if you need any pattern dental extracing forceps, We will also produce you immediately, English & American pattern



1001
Fig. 4
lower incisor



1001
Fig. 4
lower incisor



1002
Fig. 8
lower premolar



1002
Fig. 8
lower premolars



1003
Fig. 1
upper centrals & canines



1004
Fig. 2
upper laterals & canines



1005
Fig. 3
upper incisors & canines



1006
Fig. 4
upper incisors & canines



1007
Fig. 7
upper premolars



1008
Fig. 8
lower premolars



1009
Fig. 13
lower premolars




1010
Fig. 13S
lower premolars
for children.


1011
Fig. 17
upper molars, right



1012
Fig. 18
upper molars, left



1013
Fig. 18A
upper molars



1014
Fig. 18B
upper molars



1015
Fig. 19
upper third molars



1016
Fig. 20
upper third molars



1017
Fig. 21
lower molars



1018
Fig. 22
lower molars



1019
Fig. 22S
lower molars
for children



1020
Fig. 22G
lower molars



1021
Fig. 23
lower molars right



1022
Fig. 24
lower molars left



1023
Fig. 29
upper roots



1024
Fig. 29S
upper roots for children



1025
Fig. 30
upper roots



1028
Fig. 30S
upper roots
for children




1026
Fig. 31
lower roots




1027
Fig. 32
lower molars




1029
Fig. 33
lower roots



1030
Fig. 33A
lower roots




1031
Fig. 33B
lower roots




1032
Fig. 33L
lower roots



1033
Fig. 37
upper incisors
for children




1034
Fig. 33S
lower roots
for children



1035
Fig. 38
lower incisors &
canines for children



1036
Fig. 39
upper molars
for children




1037
Fig. 39R
upper molars right
for children



1040
Fig. 39L
upper molars left
for children




1038
Fig. 39A
upper molars
for children




1039
Fig. 40
lower molars
for children




1041
Fig. 41
upper roots



1042
Fig. 44
upper roots




1043
Fig. 45
lower roots




1044
Fig. 48
lower roots




1045
Fig. 49
upper roots



1046
Fig. 50
upper roots






1053
Fig. 54
for separating
upper molars




1054
Fig. 55
for separating
upper molars



1055
Fig. 56
for separating
lower molars




1056
Fig. 65R
upper molars
right



1057
Fig. 65L
upper molars left



1058
Fig. 67
upper third molars




1059
Fig. 67A
upper third molars



1060
Fig. 67 1/2R
upper third molars
right




1061
Fig. 67 1/2L
upper third molars
left



1062
Fig. 68
lower roots



1063
Fig. 68A
lower roots



1064
Fig. 73
lower molars



1065
Fig. 73R
lower molars right



1066
Fig. 73L
lower molars left




1067
Fig. 73A
lower molars



1068
Fig. 73S
lower molars



1069
Fig. 73B
lower molars for children



1072
Fig. 73K
lower molars
for children



1070
Fig. 74
lower roots




1071
Fig. 74D
lower roots




1073
Fig. 74M
lower roots



1074
Fig. 74N
lower roots



1075
Fig. 74XN
lower roots and
crowded incisors



1076
Fig. 75
lower premolars



1077
Fig. 76
upper roots




1078
Fig. 76N
upper roots




1079
Fig. 76S
upper roots



1080
Fig. 79
lower third molars



1081
Fig. 79A
lower third molars



1082
Fig. 81
lower premolars
and roots



1083
Fig. 86A
lower molars



1084
Fig. 86B
lower molars



1085
Fig. 86C
lower molars



1086
Fig. 87
lower molars




1087
Fig. 91
for separating
lower molars



1088
Fig. 99
lower molars



1089
Fig. 101
upper premolars
and roots



1090
Fig. 100
upper roots



1091
Fig. 99 1/4
lower molars



1092
Fig. 99 1/2
lower molars



1093
Fig. 102
upper molars right



1094
Fig. 103
upper molars left




1095
Fig. 105
lower incisors




1096
Fig. 106
lower incisors



1097
Fig. 103
upper canines



1098
Fig. 112
upper roots



1099
Fig. 113
upper roots



1101
Fig. 136
upper premolars
and roots




1102
Fig. 137
lower incisors
and roots



1100
Fig. 114
upper roots



1103
Fig. 143
lower roots





1110
ROUTURIER (Fig. 212)
(22 ½ L)
lower molars
and wisdom left



1111
ROUTURIER (Fig. 213)
(22 ½ R)
lower molars
and wisdom right



1112
ROUTURIER (Fig. 214)
(22 ½ L)
lower molars
and wisdom left



1113
ROUTURIER (Fig. 215)
(33 ½ R)
lower molars
and wisdom right



1114
RESECTION FORCEPS
(216/Witzel)
upper roots



1115
Fig. 5



1116
Fig. 6



1117
Fig. 9



1118
Fig. 10



1119
Fig. 11



1120
Fig. 12



1121
Fig. 14



1122
Fig. 15



1123
Fig. 25



1124
Fig. 26



1125
Fig. 27



1126
Fig. 28



1127
Fig. 34



1128
Fig. 35



1129
Fig. 42



1130
Fig. 43



1131
Fig. 44A



1132
Fig. 46



1133
Fig. 47



1134
Fig. 50A



1135
Fig. 74R



1136
Fig. 74L



1137
Fig. 75R



1138
Fig. 75L



1139
Fig. 94



1140
Fig. 95



1141
Fig. 98




1142
Fig. 105




1143
Fig. 106




1144
Fig. 2002




1145
Fig. 1
upper laterals
and canines



1146
Fig. 2
upper laterals
and canines




1147
Fig. 3
upper incisors
and canines



1148
Fig. 4
lower incisors
and canines



1149
Fig. 7
upper premolars



1150
Fig. 13
lower premolars




1151
Fig. 13S
lower premolars
for children



1152
Fig. 17
upper molars right



1153
Fig. 18
upper molars left



1154
Fig. 18A
upper molars



1155
Fig. 18B
upper molars



1156
Fig. 19
upper third molars



1157
Fig. 22
lower molars



1158
Fig. 22S
lower molars
for children



1159
Fig. 29
upper roots



1160
Fig. 29S
upper roots
for children



1161
Fig. 30
upper roots



1162
Fig. 30S
upper roots
for children



1163
Fig. 33
lower roots



1164
Fig. 33A
lower roots



1165
Fig. 33B
lower roots



1166
Fig. 33L
lower roots



1167
Fig. 33S
lower roots
for children



1168
Fig. 37
upper incisors
for children



1169
Fig. 39
upper incisors
for children




1170
Fig. 39A
upper molars
for children




1171
Fig. 39R
upper molars right
for children




1172
Fig. 39L
upper molars left
for children




1173
Fig. 41
upper roots




1174
Fig. 44
upper roots




1175
Fig. 45
lower roots



1176
Fig. 48
lower roots



1177
Fig. 51
upper roots



1178
Fig. 51A
upper roots



1179
Fig. 51S
upper roots
for children




1180
Fig. 51L
upper roots



1181
Fig. 52
upper roots



1182
Fig. 54
for separating
upper molars



1183
Fig. 65R
upper molars right



1184
Fig. 65L
upper molars left




1185
Fig. 67
upper third molars



1186
Fig. 67A
upper third molars



1187
Fig. 67 1/2 R
upper third molars
right



1188
Fig. 67 1/2 L
upper third molars
left



1189
Fig. 73
lower molars



1190
Fig. 74
lower roots




1191
Fig. 74M
lower roots



1192
Fig. 74N
lower roots



1193
Fig. 74XN
lower roots and
crowded incisors



1194
Fig. 75
lower premolars



1195
Fig. 76
upper roots



1196
Fig. 76N
upper roots




1197
Fig. 81
lower premolars
and roots



1198
Fig. 76A
lower third molars



1199
Fig. 79
lower third molars



1200
Fig. 86A
lower molars



1201
Fig. 86B
lower molars



1202
Fig. 86C
lower molars



1203
Fig. 99
lower molars



1204
Fig. 99 1/4
lower molars



1205
Fig. 136
upper premolars
and roots



1206
Fig. 143
lower roots



1207
KLEIN (Fig. 217)
upper incisors

1208
KLEIN (Fig. 218)
upper premolars



1209
KLEIN (Fig. 219)
upper roots



1210
KLEIN (Fig. 220)
upper molars



1214
ALBRECHT
Fig. 224
Deciduous

1215
ALBRECHT
Fig. 225
Deciduous




1211
KLEIN (Fig. 221)
lower incisors



1212
KLEIN (Fig. 222)
lower molars



1213
KLEIN (Fig. 223)
lower roots




1217
WITZEL
Fig. 227
upper roots


1216
WITZEL
Fig. 226
lower roots



1218
ARTHUR
Fig. 228
root fragments



1219
STIEGLITZ
Fig. 229
root fragments




1220
STIEGLITZ
Fig. 230
root fragments




1221
RALK
Fig. 231
root fragments




1222
W1-WINTER
upper centrals
and canines




1223
Fig. 1
upper incisors and roots



1224
Fig. 2
upper and lower roots



1225
Fig. 4
upper and lower premolars



1226
PHYSICK
Fig. 5
lower third molars



1227
Fig. 6
separating
forceps



1228
Fig. 7
upper roots




1229
Fig. 8
upper and lower
third molars



1231
Fig. 10H
upper molars




1232
Fig. 13
upper incisors



1230
Fig. 10S
upper molars




1233
Fig. 14
lower incisors



1234
Fig. 16
lower molars



1235
HARRIS
Fig. 15
lower molars



1237
HARRIS
Fig. 18R
upper molars right




1238
HARRIS
Fig. 18L
upper molars left



1236
Fig. 16S
lower molars
for children



1239
Fig. 17
lower molars






1246
Fig. 30S
upper teeth
for children



1247
PARMLY
Fig. 32
upper premolars
and molars



1248
PARMLY
Fig. 32A
upper premolars
and roots



1249
Fig. 34
lower anteriors
and roots




1250
Fig. 35
upper roots




1251
Fig. 36
lower anteriors
and roots



1252
Fig. 37
lower anteriors
and roots



1255
Fig. 53L
upper molars left



1253
Fig. 38
upper anteriors
and roots



1254
Fig. 53R
upper molars right



1256
Fig. 62
incisors, premolars,
deciduous teeth




1257
Fig. 65
upper incisors
and roots



1258
TOMES
Fig. 69
upper and lower roots



1259
TOMES
Fig. 85
lower canines,
premolars and molars



1260
NEVIUS
Fig. 88R
upper molars right




1261
THOMAS
Fig. 85A
lower canines,
premolars and molars



1262
NEVIUS
Fig. 88L
upper molars left



1263
KELLS
Fig. 99A
upper incisors,
canines, premolars



1264
Fig. 104
upper third molars




1266
HULL
Fig. 101
deciduous teeth



1267
Fig. 103
lower incisors,
premolars



1265
KELLS
Fig. 99C
upper incisors,
canines, premolars



1268
CRYER
Fig. 150
lower incisors,
premolars, roots



1269
CRYER
Fig. 150A
upper incisors,
premolars, roots



1271
CRYER
Fig. 151
lower incisors,
premolars, roots



1272
CRYER
Fig. 151A
lower incisors,
premolars, roots




1270
CRYER
Fig. 150S
upper teeth for children




1273
CRYER
Fig. 151S
lower teeth
for children



1274
Fig. 203
lower incisors
and roots



1275
Fig. 210S
upper third molars



1278
Fig. 213
upper centrals and canines




1276
Fig. 201
upper incisors,
canines, premolars



1277
Fig. 210H
upper third molars



1279
Fig. 217
lower molars




1280
Fig. 222
lower third molars



1282
Fig. 286
Swing lower molars




1283
Fig. 288
lower canines
and premolars




1281
Fig. 286
upper incisors and roots



1284
COHEN
Fig. 310
universal forceps




1285
NEVIUS
Fig. 320/10S
upper canines



1286
HENAHAN
Fig. 315/1A upper incisors




1287
NEVIUS
Fig. 325/10H
upper canines




1288
STANDARD
Fig. 330/4
upper third molars




1289
WOODWARD
Fig. 335/FH
lower molars



1290
Fig. 340/89
upper molars right



1292
Fig. 345/90
upper molars left



1291
WOODWARD
Fig. 335/FS
lower molars



1293
Fig. 350/23S
lower molars



1294
Fig. 355/23SK
lower molars



1295
Forcep Fig. 301



1295-A
Forcep Fig. 300



1296
Fig. 360/17SK
lower molars



1297
Fig. 365/150SK
upper anteriors
and roots



1298
Fig. 370/151SK
lower anteriors
and roots



1299
Molt 9



1300
Allen 9A



1301
Prichard



1302
Seldin



1303
Goldman- Fox 14



1304
West
Dell W2

Heidbrink Root - Tip Picks

Root Tip



1305
DEL 78



1306
DEL 79



1307
DEL 80



1308
DEL H1



1309
DEL H2



1310
DEL H3



1311
DEL H2-3



1312
POTT
6X



1313
POTT
7X



1314
BUSER



1315
OHI 10



1316
HOWARD 9H



1317
FREER



1318
FREER 15



1319
FREER 16

lower jaw, left side

lower jaw, right side




1320
BEIN
Fig. 4



1331



1332



1321
BEIN
1.0 mm

1322
BEIN
1.5 mm

1323
BEIN
2.0 mm

1324
BEIN
2.5 mm

1325
BEIN
3.0 mm

1326
BEIN
3.5 mm

1327
BEIN
4.0 mm

1328
BEIN
4.5 mm

1329
BEIN
5.0 mm

1330
BEIN
5.5 mm



1333
BEIN



1334
BEIN
Fig. 1




1335
BEIN
Fig. 2



1336
BEIN
Fig. 0



1337
BEIN
Fig. 3



1338
BEIN
Fig. 4



1339
BEIN
Fig. 5



1340
BEIN
Fig. 6



1341
CRYER
Fig. 3



1342
CRYER
Fig. 4



1343
CRYER
Fig. 12



1344
CRYER
Fig. 27



1345
CRYER
Fig. 28



1346
CRYER
Fig. 27S



1347
CRYER
Fig. 28S



1348
SELDIN
Fig. 2



1349
SELDIN
Fig. 3



1350
SELDIN
Fig. 4L



1351
SELDIN
Fig. 4R



1352
SELDIN
Fig. 1L



1353
SELDIN
Fig. 1R



1354




1355




1356




1357



1358



1359



1360

Root Screws




1362
Morrison
2.0 mm



1363
Morrison
2.5 mm



1364
Morrison
3.0 mm



1361
Morrison
3.5 mm



**1365
IS-APICAL**



**1366
APICAL
2 mm
Fig. 1**



**1367
APICAL
3 mm
Fig. 2**



**1368
APICAL
4 mm
Fig. 3**



**1369
APICAL
4.5 mm
Fig. 4**



**1370
APICAL
2 mm
Fig. 5**



**1371
APICAL
3 mm
Fig. 6**



**1372
APICAL
4 mm
Fig. 7**



**1373
APICAL
4.5 mm
Fig. 8**



**1374
APICAL
Fig. 9**



**1375
WINTER
Fig. 14 L**



**1376
WINTER
Fig. 14 R**



**1377
HEIDBRINK
Fig. 1**



**1378
HEIDBRINK
Fig. 2**



**1379
FLOHR
Fig. 1**



**1380
FLOHR
Fig. 2**



**1381
APICAL
Fig. 1**



**1382
APICAL
Fig. 2**



1383
IS-CRYER



1384
CRYER
Fig. 1



1385
CRYER
Fig. 2



1386
IS-COUPLAND



1387
COUPLAND
Fig. 1



1388
COUPLAND
Fig. 2



1389
COUPLAND
Fig. 3



1390
Fig. 1



1391
Fig. 2



1392
Fig. 3



1393
FLOHR
Fig. 1



1394
FLOHR
Fig. 2

1395
FLOHR
Fig. 3



**1396
KOPP
Fig. 1**



**1397
KOPP
Fig. 2**



**1403
KOPP
Fig. 1**



**1404
KOPP
Fig. 2**



**1396
KOPP**



**1398
KOPP
Fig. 3**



**1399
KOPP
Fig. 4**



**1402
KOPP**



**1405
KOPP
Fig. 3**



**1406
KOPP
Fig. 4**



**1400
KOPP
Fig. 5**



**1401
KOPP
Fig. 6**



1407
HYLIN




1407
HYLIN
Fig. 1




1408
HYLIN
Fig. 2




1409
BEIN




1409
BEIN
Fig. 1



1411
BEIN
Fig. 3



1410
BEIN
Fig. 2



1412
BEIN
Fig. 4




1413
WARWICK-JAMES
Fig. 1



1413
Warwick-James
Fig. 2



1414
Warwick-James
Fig. 3



1415
London-Hospital
Fig. 1



1415
London-Hospital
Fig. 2



1416
London-Hospital
Fig. 3



1417
OTT



1418
CHOMPRET
Fig. 1



1419
CHOMPRET
Fig. 2



1420
CHOMPRET
Fig. 3




1421
CHOMPRET
Fig. 4



1423
Fig. 1



1422
Fig. 2



1424
Fig. 1



1425
Fig. 2



1426
APICAL



1427
IS-Fig. 1



1428
Fig. 2



1429
Fig. 3



1430
HEIDBRINK
Fig. 1



1431
HEIDBRINK
Fig. 2



1432
HEIDBRINK
Fig. 3



1438
MEAD
MD1



1433
LECLUSE



1434
MEAD
MD-5 A



1435
MEAD
MD-6 A



1436
MEAD
MD-7 A



1437
MEAD
MD-8



1439
WINTER
Fig. 1L




1440
WINTER
Fig. 1R



1441
WINTER
Fig. 2L



1442
WINTER
Fig. 2R



1439
WINTER



1443
WINTER
Fig. 11L



1444
WINTER
Fig. 11R



1445
WINTER
Fig. 12L



1446
WINTER
Fig. 12R



1447
WINTER
Fig. 14L



1448
WINTER
Fig. 14R



1449
BARRY
Fig. 1



1450
BARRY
Fig. 2



1451
BARRY
Fig. 3



1452
BARRY
Fig. 4



1449
BARRY



1453
BARRY
Fig. 11



1454
BARRY
Fig. 12



1455
BARRY
Fig. 17



1456
BARRY
Fig. 18



1458
CRYER
Fig. 26



1459
SELDIN
Fig. 27



1460
CRYER
Fig. 28




1461
CRYER
Fig. 29



1462
HEIDBRINK
Fig. 22



1463
HEIDBRINK
Fig. 23



1457



1465
KOPP
Fig. 31



1466
KOPP
Fig. 32



1467
KOPP
Fig. 33




1468
KOPP
Fig. 34



1469
KOPP
Fig. 35



1470
KOPP
Fig. 44



1464



1472
HEIDBRINK
Fig. 10



1473
WINTER
Fig. 11



1474
WINTER
Fig. 12



1475
BARRY
Fig. 17



1476
BARRY
Fig. 18



1477
BARRY
Fig. 20

Root Elevators according to Dr. Sader,
with interchangeable working tips



1471



1478




1479




1480

1481
Root Elevator, Plastic Handle, Set of 10 Pcs



1482
Root Elevator, Plastic Handle, Set of 5 Pcs





Manufacturer & Supplier Premium Quality Surgical, Dental,
Podiatry, Pet Grooming and Beauty Care Instruments.

USA Address:

 5900 Balcones Drive # 16244
Austin TX, USA 78731.
 +1 (917) 773 3694
 **uk:** +44 7871 836765
 info@neoaindustry.com
neoaindustry4@gmail.com
 www.neoaindustry.com

Main Head Office:

 Roars Road New Miana Pura
West, Sialkot 51310 Pakistan.
 +92 322 5587688
 +92 322 5587688
 info@neoaindustry.com
neoaindustry4@gmail.com
 www.neoaindustry.com

